



iM8 Series

Patient Monitor



iM8 Series
Patient Monitor




iM8 Series

iM8/iM8A Reliable Patient Monitor



Rolling Stand Wall Mount

The iM8 series put reliable monitoring technologies at your finger tips. With their accurate readings and optimized performance, the iM8 series can meet basic monitoring needs, bringing easy access to high-quality healthcare for everyone.



iM8 (12.1")

iM8A (10.4")

- High resolution color TFT-LCD screen
- Pacemaker detection
- Defibrillation and electrosurgical interference protection
- Pitch tone (Pulse-tone modulation)
- Unique iSEAP algorithm which is specially optimized for arrhythmia patients
- Dual-mode anti-interference pulse oximetry
- SP10 Validated NIBP
- Built-in rechargeable Li-ion battery
- Nurse call
- VGA output (optional)
- Suitable for adult, pediatric and neonatal patients


Standard Parameters: 3/5-lead ECG, RESP, EDAN SpO₂, NIBP, PR, 2-TEMP

Optional: 2-IBP, Respironics CO₂, EDAN G2 CO₂, Thermal Recorder, WLAN Accessory Kit

500 NIBP Measurement Review	96 h Trend Review	60 Alarm Review	720 s Frozen Waveform
---------------------------------------	-----------------------------	---------------------------	---------------------------------

Seamless Connectivity

The Wi-Fi connectivity of the iM8 series enables healthcare providers to monitor their patients' health status from almost anywhere. Connecting it with EDAN MFM-CMS software, you may log on from any PC/Tablet to and check the status of your patients.



End tidal Carbon Dioxide for Intubated/Non-Intubated Patients



Respironics CO₂ (Mainstream/Sidestream)

- Plug & play module design
- Dehumidification tube instead of water trap
- Low sampling rate of 50ml/min suitable for all types of patients

EDAN G2 CO₂ (Sidestream). NO FDA Respironics CO₂ FDA APPROVED.

- Superior water trap design for accurate monitoring
- iCARB™ algorithm with intelligent CO₂ pseudo wave identification technology
- Multiple sampling accessories as options for adult, child and neonate patients

Display Modes



3/5-lead ECG Analysis

OxyCRG for Neonate

Large Font

Trend Screen